


About Gator Bio

Gator Bio develops, manufactures, and markets life science analytical technologies including Gator[®] systems based on the next-gen Biolayer Interferometry. The company was founded by the industry veterans Dr. Hong Tan and Mr. Bob Zuk. Previously, Dr. Hong Tan founded ForteBio[®] and led the invention of Octet[®] BLI technology.

Gator Bio together with its sister companies have more than 600 employees worldwide and sell both diagnostics and research-use-only products. The company is ISO13485 certified. Gator[®] systems have been adopted by scientists and researchers in North America, Asia Pacific, Europe, and Middle East.

The investors of the company include Legend Capital, Matrix Partners, Maison Capital, Qiming Venture, HillHouse, Sequoia Capital, Kaiser Permanente, and Sinovation etc.

Gator Bio, Inc.
2455 Faber Place
Palo Alto, CA 94303, USA
+1 855 208 0743
info@gatorbio.com
© 2022 Gator Bio, Inc. All rights reserved. Gator is a registered trademark of Gator Bio, Inc.



General Specifications

	Gator [®] Prime	Gator [®] Plus	Gator [®] Pro
Detection Technology	Next-gen Biolayer Interferometry		
Simultaneous Reads	8	8	8, 16, 24, and 32
Spectrometers	8	8	32
Acquisition Rate	2, 5, and 10 Hz		
Temperature Control	Ambient plus 4°C to 40°C		
Dimension (HxWxD) and Weight	47 x 31 x 67 cm, 35 kg	68 x 44 x 73 cm, 55 kg	84 x 114 x 77 cm, 220 kg
Automation Compatible	No	No	Yes

Performance Specifications

	Gator [®] Prime	Gator [®] Plus	Gator [®] Pro
Sample Types	Proteins, antibodies, peptides, nucleic acids, liposomes, viruses, small molecules		
Maximum Sample Capacity	168	456	1152
Types of Analysis	Yes/no binding, quantitation, kinetics, affinity, off-rate ranking, epitope binning		
Minimum Sample Volume	100 µL	40 µL	40 µL
Baseline Noise (RMS)	≤ 4 pm		
Baseline Drift	≤ 0.12 nm/hour	≤ 0.1 nm/hour	≤ 0.1 nm/hour
Associate Rate k _{on}	10 ¹ to 10 ⁷ M ⁻¹ s ⁻¹		
Dissociation Rate k _{off}	10 ⁻⁶ to 10 ⁻¹ s ⁻¹		
Affinity Constant K _D	10 pM to 1 mM		
Binning Capacity	12x12	16x16	32x32


gator



Gator[®] Label-Free Analysis Systems

The Next-Gen Biolayer Interferometry

GatorBio.com



Biolayer Interferometry (BLI)

A Powerful Tool for Discovery, Development, and Manufacturing


A User-Friendly Label-Free Technology

A Full Suite of Applications

Gator® systems are label-free analysis instruments based on next-gen biolayer interferometry (BLI) technology. BLI detects biomolecular interactions by immersing biosensing probes in samples.

Gator® probes are micro glass rods with the distal ends coated with proprietary optical layers and surface chemistries.

The association or disassociation of biomolecules causes a phase-shift of the optical interference pattern generated from a probe’s sensing surface. Continuous measurements of the phase-shift yield binding curves.



The sensorgram shows the real-time association and disassociation curves for a binding kinetics experiment using a Gator® system.

The ease of use, versatility, flexibility, and throughput of Gator® systems have enabled many applications in therapeutic development, manufacturing, and life science research.

The next-gen BLI demonstrates higher sensitivity and more robust performance than the other commercial BLI products. It also supports wider range of applications, from drug discovery to therapeutics manufacturing.

- Biotherapeutics**
 - Antibody titer measurements
 - Kinetics analysis
 - Epitope binning
 - Process development
 - Manufacturing QC
 - Pharmacokinetics
- Gene Therapy**
 - AAV quantitation & kinetics
 - Receptor interaction
 - Gene expression
 - Neutralizing/ Total Antibody Detection

- Drug Discovery & Development**
 - Protein - small molecule interaction
 - Peptide binding analysis

- Life Science Research**
 - Protein - protein interaction
 - Receptor - ligand binding
 - Assay development and optimization

Gator® Systems consist of instruments, probes, and integrated data acquisition and analysis software pack-age.

- Simple and fast assay setup
- Automated quantitation
- Quantitation, kinetics, and regeneration in one run
- Kinetics and affinity analysis
- Real-time binding curves
- Epitope binning
- Assay template generation
- Report generation

Gator® Software for GMP and GLP

Gator® Part11 Software enables users in GMP or GLP environments to comply with FDA 21 CFR Part 11 regulations. All data acquired with the Part11 Software is time-stamped and traceable. Features such as account management, enhanced audit trails, and recorded user sessions are in compliance with FDA guidance.

Gator® Probe	Function	Application	Dynamic Range	Regeneration
ProA	IgG titer	Q	0.02 – 2000 µg/mL	Yes
ProG	IgG titer	Q	0.02 – 2000 µg/mL	Yes
ProL	Kappa light chain titer	Q	0.02 – 2000 µg/mL	Yes
SA	Biotinylated and Avi-tagged molecules	K/EP	Protein dependent	No
SA XT	Biotinylated proteins and large molecules	K	Protein dependent	No
Flex SA Kit	Reusable SA kit	K	Protein dependent	Yes
SMAP	Measurement of small molecules, peptides (<150 Da)	K	Protein dependent	No
HFC	Human IgG characterization	Q/K/QKR/EP	0.05 – 300 µg/mL	Yes
HFCII	Advanced human IgG characterization	Q/K/QKR/EP	0.3 – 6000 µg/mL	Yes
MFC	Mouse IgG characterization	Q/K/QKR/EP	0.02 – 2000 µg/mL	Yes
Anti-FAB	F(ab), F(ab')2	Q/K/QKR/EP	0.3 – 3000 µg/mL	Yes
APS	Direct adsorption	K	Protein dependent	No
AR	Amine coupling immobilization	K/EP	Protein dependent	No
His	His-tagged proteins	Q/K/QKR/EP	Protein dependent	Yes
Ni-NTA Kit	His-tagged proteins through Ni-NTA	Q/K/QKR/EP	0.25 – 1000 µg/mL	Yes
AAVX	Direct binding titer (AAV1-10)	Q/K	1x10 ⁹ – 1x10 ¹³ vp/mL	Yes
AAV9	Direct binding titer (AAV9)	Q/K	3x10 ⁹ – 1x10 ¹³ vp/mL	Yes
HS AAV Kit	High sensitivity titer (AAV1-8, 10)	Q	1x10 ⁷ – 5x10 ¹⁰ vp/mL	No
HS AAV9 Kit	High sensitivity titer (AAV9)	Q	1x10 ⁷ – 1x10 ⁹ vp/mL	No
AAV Ratio Kit	Empty vs Full Ratio Determination	Ratio	0-100% full	No
Anti-Rabbit	Rabbit-IgG	Q/K	1 – 500 µg/mL	Yes
Anti-FLAG	FLAG-tagged proteins	K	Protein dependent	No
SARS-RBD	RBD antibodies and receptors	Q	Protein dependent	Yes