

# Evaluating the Versatility of Gator® Next Generation BLI Platform for Biotherapeutic Development and AAV Analytics

Harsha Agarwal, Johnny Zhang, Wai Choi, Lu Peng, Indrani Chakraborty, Vivienne Lee  
Gator Bio, Palo Alto, CA

## Introduction


### Bi-layer Interferometry (BLI)

- Label-free technology based on reflection of light on the surface of a biosensor tip
- The shift in interference pattern plotted against time when a molecule is bound
- The change in pattern proportional to the number of biomolecules bound
- Gator® next-generation BLI is a versatile real-time analysis platform
- Minimal hands-on time
- Wide applications ranging from protein-protein interactions, therapeutics development and viral vector analysis
- Tolerant to different buffers, cell media, crude lysates, serum and plasma

Here, we present data from some unique applications such as,

- Epitope binning
- Biosimilar kinetics
- LNP quantitation
- Nanobody screening
- Small molecule interactions
- AAV analytics


Gator® Bio's broad portfolio of biosensors can support at multiple stages of therapeutic development and in gene therapy.



## Small molecule interactions


### Gator® SMAP Probes

- Detect small molecules up to 150 Da with a binding partner
- Streptavidin based proprietary surface chemistry for high-capacity immobilization of biotinylated binding partner
- Determination of the kinetic parameters ( $k_{on}$ ,  $k_{off}$  and  $K_D$ ) of the small molecule with the immobilized binding partner
- Enhanced signals vs traditional BLI platforms





### Peptide binding using Gator® Streptavidin (SA) XT Probes

- The biosensors detect biotinylated oligos, peptides and proteins above 1 kDa
- Unique optical layer with novel proprietary chemistry enhances the signal 5-3x than the traditional BLI platform
- Large biomolecule up to 2 MDa can be detected without inversion of signal
- Higher signal allows for lower loading of ligand and analyte, hence conserving precious sample
- Accurate determination of the kinetic parameters ( $k_{on}$ ,  $k_{off}$  and  $K_D$ ) from small peptides to large biomolecules




## Epitope Binning

- High-throughput 32 x 32 epitope binning assay utilizing Gator® Pro instrument in less than 8 hours
- Accurate and automated tandem or traditional sandwich format
- Easy data visualization and interpretation
- Broad range of biosensors for tandem and sandwich format




## Biosimilar Kinetics

- Gator® Human Fc (HFC) Receptor II biosensor detects Fc region of all four IgG isotype making it suitable for biosimilar screening
- The biosensor can be regenerated up to 20 times without loss in signal, thus making them cost effective
- No cross-reactivity to Human Fab region and other species antibody, making them very specific



## Lipid Nanoparticle Quantitation


- Anti-PEG probes enables the detection and quantitation of LNPs
- Using unique optical layer, no inverted binding signal from LNP binding
- Serum proteins can be immobilized onto the probes to study the interaction with LNPs



## Nanobody Screening

- One-of-a-kind nanobody screening biosensor (Anti-VHH).
- Fast and accurate twin-strep-tagged nanobody screening using Strep-Tactin® XT biosensors.
- $K_D$  comparison studies of nanobody screening using both the biosensors shows comparable results.

### $K_D$ Comparison with Anti-Strep Tag Biosensors



## Conclusion

- Fast biomolecule characterization using Gator® Pro Instrument
- Automated, accurate and fast epitope binning
- Precise and efficient biosimilar kinetics
- Accurate and easy kinetic platform for LNPs
- Easy and specific nanobody screening
- Enhanced small molecule and protein kinetic interactions
- Accurate kinetic parameters for small peptides to large protein using Gator® SA XT biosensors
- Total AAV solutions:
  - Precise AAV titer from upstream to downstream samples
  - Accurate determination of Empty vs full ratio
  - Easy and fast detection of Genome titer

