



Turn-Key Automation Solution with Gator BLI

For use with Gator® Pivot and Gator® Pro

Introducing, the Gator Automation Solution. Engineered for seamless integration with a variety of lab equipment, the Gator Automation Solution offers a simple way to configure Gator BLI instruments to fit your workflow — enabling automation across multiple levels of your lab operations.

Integrate with:


Gator Automation Solution



Liquid Handler



Incubator



**Biosensor & Sample
Plate Hotel**



BLI Instrument





Gator Automation Solution


Flexible, Easy, Powerful



Compatible with 96-well and
384-well SBS plates



Load and unload sample plates
and biosensor plates



Uncover and recover plates to
avoid evaporation

Get in touch with us
We're always here to help

1-855-208-0743

info@gatorbio.com

gatorbio.com

